

**Germ Thief** is a non-woven, disposable wipe, pre-saturated with a quaternary disinfectant cleaner. **Germ Thief** is intended for use in:

• Hospitals and healthcare settings: ambulances, anesthesia, CAT labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctor's office(s), Emergency Rooms, E.R., Home health care, Hospices, ICU, Laboratory, Laundry rooms, Neonatal Intensive Care Units (NICU), Newborn nursery, Nursing homes, Operating rooms, Pediatrics, Pediatric Intensive Care Units (PICU), Physical therapy, Physician's office(s), Ophthalmic offices, Radiology, Recovery Room, Rehabilitation, Respiratory therapy, Surgery, Transport vehicles, X-Ray

• Critical Care Areas: CCU, Emergency Rooms, E.R., I.C.U. Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU), Surgery, Ophthalmic offices

• Schools, health clubs, gyms, athletic facilities, airports, transportation terminals **Germ Thief** will clean and disinfect hard, non-porous, surfaces of: bed railings, cabinets, carts, chairs, counters, exam tables, gurneys, isolettes, IV poles, stethoscopes, stretchers, tables, telephones, toilet seats, toilets, bathtubs, bathroom fixtures, shower stalls, sinks, tiles, vanity tops, hospital beds, traction devices, MRI, CAT, scales, paddles, wheelchairs, ultrasound transducers and probes. **Germ Thief** will clean and disinfect hard, non-porous surfaces of: gym equipment, exercise machines, wrestling mats, gymnastic equipment.

#### For Ultrasound Transducers and Probes:

Use only on external surfaces of ultrasound transducers and probes. With a clean washcloth, clean gel from external surface of the ultrasound transducer or probe to be disinfected. Then thoroughly wet probe with wipe. Treated surface must remain wet for 10 minutes. Use additional wipes if needed to assure continuous 10 minute wet contact time. Let air dry.

**KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings (Hospitals, Nursing Homes) or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS).

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS. PERSONAL PROTECTION:** When handling items soiled with blood or body fluids use disposable latex gloves, gowns, masks, and eye coverings.

**CLEANING PROCEDURES:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of this product. **Germ Thief** wipes can be used for this purpose. **CONTACT TIME:** Thoroughly wet surface with a new **Germ Thief** wipe. Allow surface to air dry. Although efficacy at a 2 minute contact time has been shown to be adequate against HIV-1, use the contact time for the disinfection of other organisms listed on this label.

**DISPOSAL OF INFECTIOUS MATERIALS:** Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal. Not for personal cleansing • This is not a baby wipe • Do not flush in toilet.



**Cleans, Disinfects & Deodorizes**

**Kills 99.9% of Bacteria**

**Bleach Free • Alcohol Free**

**Tuberculocidal • Virucidal**

**Pseudomonocidal • Bacterial**

**Kills Pandemic 2009 H1N1 Influenza A Virus**

#### Active Ingredients

n-Alkyl (60% C<sub>14</sub> 30% C<sub>16</sub> 5% C<sub>18</sub>) dimethyl benzyl ammonium chlorides.....0.105%

n-Alkyl (68% C<sub>12</sub> 32% C<sub>14</sub>) dimethyl ethylbenzyl ammonium chlorides.....0.105%

INERT INGREDIENTS.....99.790%

TOTAL.....100.000%

EPA REG. NO. 1839-83-89808

EPA EST. NO. 8155-OH-1

KEEP OUT OF REACH OF CHILDREN  
CAUTION

See Side (back) Panels for additional Precautionary Statements

**Net Contents 195 6x12 Pre Moistened Disinfectant Wipes**

P.I.N. D-116

US Patent No. 5,444,094

ENCORE INDUSTRIAL PRODUCTS, LLC

281-944-4777 • 1-877-703-3044

3502 Bacor Rd. Houston, TX 77084

www.encoreindustrial.com

#### DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body.

**CLEANING:** Use one **Germ Thief** wipe to completely remove surface of all gross debris from surfaces to be treated.

**DISINFECTING:** To disinfect hard, non-porous surfaces, thoroughly wet the surface to be treated. Treated surface must remain visibly wet for disinfection:

• To disinfect against *Escherichia coli* (*E. coli*), *Escherichia coli* O157:H7, *Pseudomonas aeruginosa*, *Salmonella* (*choleraesuis*) *enterica* and Vancomycin Resistant *Enterococcus faecalis* (VRE), allow the treated surface to remain wet for 10 minutes.

• To disinfect against *Mycobacterium bovis* BCG (Tuberculosis), allow the treated surface to remain wet for 10 minutes.

• To disinfect against *Staphylococcus aureus* and Methicillin Resistant *Staphylococcus aureus* (MRSA), allow the treated surface to remain wet for 10 minutes.

• To disinfect against Vancomycin intermediate Resistant *Staphylococcus aureus* (VISA), allow the treated surface to remain wet for 10 minutes. A potable water rinse is required for surfaces in direct contact with food. PUSH FINGER THROUGH HOLE. Pull through about one inch. Replace cover. Pull out first wipe and snap off at a 90 degree angle. Remaining wipes feed automatically ready for next use. When through using, keep small center cap closed to prevent moisture loss. Do not flush.

**DISPENSER:** To start feed. Remove large cover and discard seal from container. From center of wipe roll, pull up wipe corner, twist into a point and thread it through the hole in the container cover

#### PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals

**CAUTION:** Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum or using tobacco.

#### FIRST AID

Have the product container or label with you when calling a poison control center or doctor, or going for treatment. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. **IF IN EYES:** Hold eyes open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.

#### STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

**PESTICIDE STORAGE:** Store in an area inaccessible to children.

**CONTAINER STORAGE/DISPOSAL:** Wipes: Dispose of wipe in trash after use. Do not flush. Canister: Tightly close lid between uses to retain moisture. Nonrefillable container. Do not reuse or refill this container. Offer empty container for recycling if available or discard in trash.



Consult Directions